Red Light Therapy Panel Max – GymBeam
Red Light Therapy Panel Max – GymBeam
Pannello per Terapia a Luce Rossa Max – GymBeam
Panel do terapii światłem czerwonym Max – GymBeam
Red Light Therapy Panel Max – GymBeam
Red Light Therapy Panel Max – GymBeam
Vörösfény-terápiás Panel Max – GymBeam
Red Light Therapy Panel Max – GymBeam
Red Light Therapy Panel Max – GymBeam
Панел за терапия с червена светлина Max – GymBeam
Панель для червоної світлотерапії – GymBeam
Red Light Therapy Panel Max – GymBeam
Red Light Therapy Panel Max – GymBeam
Red Light Therapy Panel Max – GymBeam

GymBeam

User Manual

Technical Specifications

Features	
LED Quantity	432 Pcs (5W Dual Chip)
Input Voltage	100–240V AC
Dimensions	92 x 42 x 6.5 cm
Spectrum	630, 650, 660, 670, 810, 830, 850, 1060 nm
Irradiance	165 mW/cm^2 (from a distance of 15 cm)
Timer	1–30 Mins
Frequency	0–10,000 Hz
Brightness	0–100%
Frame Material	Steel


Safety Instructions

Please read this manual carefully before installation and keep it for future reference.

1. **Children:** This device is not intended for use by children. Always consult a physician before use on children.
2. **Medical Conditions:** Consult your doctor before use if you are pregnant, have a pacemaker, or other health concerns.
3. **Photosensitivity:** If you take medications that increase light sensitivity or have a photosensitive skin condition, consult a doctor first.
4. **Eye Protection:** Use of the provided protective goggles is required during therapy. Do not look directly into the LEDs.
5. **Exposure:** The treated skin area must be completely bare (without clothing) during use.
6. **Adverse Reactions:** If itching or redness occurs within 48 hours, stop using the device and consult a physician.
7. **Damage:** Do not use the device if any part is damaged (broken timer, missing LEDs, etc.).
8. **Repairs:** Do not attempt to repair or modify the device yourself. Repairs must be performed by a qualified professional.
9. **Environment:** Indoor use only. Do not expose to moisture, dust, or water sources.
10. **Purpose:** This device is not a substitute for professional medical care.

Product Description

1. **Near-Infrared Light (NIR):** LEDs emitting invisible spectrum light.
2. **Red Light (R):** LEDs emitting visible red light.
3. **LED Chips:** 432 high-power chips.
4. **Touch Screen:** Interactive display for settings and control
5. **Hanging System:** Pulley and cables for door mounting.

Package Contents

- 1x Red Light Therapy Panel Max
- 1x Power Cable
- 1x User Manual
- 1x Protective Goggles (Glasses)
- 1x Remote Control
- 1x Hanging Set (Cables, Pulley, Door Hooks)
- 2x Tandem Screws

Operating Instructions

The device is operated either with the remote control or directly on the display.

Main Panel

1. Time and mode display.
2. Scene: Press to select a pre-set setting that best suits your needs (Fat Reduction / Skin / Muscles / Joint Pain / Universal).
3. Start/Pause: Press to turn the light on or off.
4. Settings: Click to customize the timer, mode, brightness, frequency, group device linking, alarm settings, date and time, or other settings.

Setting Options

1. Time Setting

Supports time settings from 0 to 30 minutes. The time can be set in two ways:
- Slide the bar with your finger around the circle to adjust the time
- Tap the centre of the circle to go to the next page, where you can enter the desired time and tap the tick symbol to save it.

2. Mode Selection

There are 3 modes available: RED (red light), NIR (near-infrared light), and RED+NIR (combination).
- RED: Red light, operates at wavelengths of 630 nm and 660 nm.
- NIR: Near-infrared light, operates at wavelengths of 810 nm, 830 nm, and 850 nm. Please note that **NIR light is invisible to the human eye.**

3. Brightness Setting

Supports 0-100% brightness. The brightness can be set in two ways:
- Tap the "+" or "–" buttons on the left to adjust the brightness.
- Move the slider icon on the right side of the screen to adjust the brightness level.

4. Frequency

Supports a frequency from 0 to 10,000 Hz. The frequency can be set in two ways:
- Tap the "+" or "–" buttons to adjust frequency.
- Tap the "Null" button to open a numeric keypad. Enter the desired frequency number and confirm.

Maintenance

1. Before cleaning, ensure the device is disconnected from the power supply.
2. Wipe the panel (except the lenses) with a damp microfibre cloth.


WARNING: Do not use chemical cleaning agents. Do not spray or apply any liquid directly onto the device, especially near the ventilation openings or LED diodes.


Storage

The device and its accessories should be stored in a clean, dry place.
- Temperature: 0 °C to 40 °C (32 °F to 104 °F)
- Humidity: less than 50%

INSTALLATION (OPTIONAL)


- **STEP 1:** Turn the face of your panel towards a soft, scratch-free surface. Screw the threaded posts into the top of the device. Next, screw the braided cables securely onto the threaded posts.
- **STEP 2:** Hang the door hooks on any standard-sized door.


GymBeam


- **STEP 3:** Hang your red light panel on the door by connecting the braided cables to the snap link of the pulley system.


- **STEP 4:** To raise the height of your red light, support the device with one hand and pull down the loose cable with your opposite hand. Ensure you adequately support the device while raising it.

- **STEP 5:** To lower your red light, push up on the lever of the height-adjusting pulley. Support the device with your opposite hand to ensure it is adequately supported while lowering.

- **STEP 6:** Locate the male power socket on the backside of the device. Insert the power cord into the socket and plug the other end into a wall outlet. Turn the power switch to the “on” position.

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